



Pfizer Vaccine Box Freezer Racks

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

Smart storage solution designed for Pfizer's Covid-19 Vaccine Boxes An ideal tool to keep vaccines organized and ready for retrieval



UFD-252-T2
Stores (10) Pfizer
Covid-19 Vaccine Boxes



UFD-352-T6
Stores (15) Pfizer
Covid-19 Vaccine Boxes

All Crystal freezer racks come with our unconditional lifetime warranty



Crystal Freezer Racks are made from pharmaceutical grade stainless steel and feature double welds. All edges are smoothed for easy handling.

Part #	Dimensions (D x W x H) (in)	Dimensions (D x W x H) (mm)	Capacity
UFD-252-T2	18.42 x 9.29 x 9.45 (in)	468 x 236 x 240 mm	10 Boxes
UFD-352-T6	27.48 x 9.29 x 9.45 (in)	698 x 236 x 240 mm	15 Boxes