

Certificate of Assurance

AltemisLab certifies that this product meets the following criteria.

Product: AlteTube™ 2D SBS format barcoded tubes

Component material: Medical Grade Virgin Polypropylene

Treatment: Vortexing and Centrifugation

Details:

Centrifugation, max. RCF

5800 x g: swinging-bucket rotor

26000 x g: fixed-angle rotor

PLEASE NOTE: For individual tubes or racked tubes, tubes must be supported from underneath to prevent damage to the rack.

Centrifugation to be done at ambient temperatures.

A handwritten signature in black ink, appearing to read 'L Frost'.

Lynsey Frost
Director

AltemisLab Ltd

Vanilla House, Copley Hill Business Park,
Babraham, Cambridge, CB22 3GN
United Kingdom
+44 (0) 1223 344900
Company Number: 11188771

